

# ORTEC®

## ANGLE 4

Advanced Gamma Spectroscopy  
Efficiency Calibration Software




“Compatible, Efficient, and Defendable Calibrations  
for Gamma Spectroscopy Applications.”

**AMETEK®**  
ADVANCED MEASUREMENT TECHNOLOGY

# ANGLE 4

ANGLE is an advanced efficiency calculation application for High Purity Germanium and Sodium Iodide detectors based on the concept of Efficiency Transfer. This method combines the measured efficiency of a known reference configuration and solid angle models to derive the efficiency for different containers, sample materials, and sample positions. This semi-empirical approach is more accurate than pure mathematical models due to large errors that can be imposed by detector characteristics that are not precisely known – such as crystal defects, contact thickness, and dead layers – as these errors cancel out in the reference efficiency measurement. And, since the Reference Efficiency can be determined from any standard source, there is no need for complex and costly factory characterization of the detector.



## WHY ANGLE?

### Compatibility

- All 32-bit and 64-bit versions of Windows from Windows 95 to Windows 10
- Multiple Language: English, French, Spanish, Russian, Chinese, Japanese – and growing!
- High Purity Germanium and Sodium Iodide Detector Types
- Modeling for most common laboratory measurement containers

### Process Efficiency

- No Factory Detector Characterization Necessary
- Import and Export ORTEC and Canberra file formats
- Command line scripting and XML Data files for automation and application integration
- Rapid modeling using Container, Geometry, and Source Matrix configurations

### Defendable Results

- Comprehensive Reporting of Efficiency Calculation Model
- Actual and Relative Efficiency Method provides Calibration Standard Traceability
- Graphical Display of model
- High Accuracy with Extensive Comparison Testing

## New Features in ANGLE 4!

- New!** Modeling for Cylinder and Well Sodium Iodide Detectors in addition to Germanium
- New!** Multi-Language Support which is easily extended to additional languages
- New!** Comprehensive and Summary Modeling Reports for validation and record retention
- New!** Graphical Display of Model and Calculated vs Reference Efficiency curves
- New!** Geometry Correction Files to maintain traceability and calibration uncertainty in GammaVision analysis
- New!** Command Line Operation and XML File Format for automation and integration with other applications
- New!** Import and Export any parameter from Efficiency Calculation Parameter or Results files
- New!** Discrete Reference Efficiency eliminates the curve fit error at specific energies of interest
- New!** Enhanced User Interface for a more intuitive user experience
- New!** Demo Detectors to evaluate and perform training with no registration required

# ANGLE 4

## ANGLE Overview

### Detector, Container, Matrix, and Position Models

#### Detector

**GEM40-83-ICS\_55-P42871A**  
GEM40-83-ICS\_55-P42871A

#### Container

**- Example Marinelli**  
The example of a Marinelli container

**1L Marinelli (132G)**  
1L Marinelli (132G) filled to 105 mm Water

#### Source

**Source height: 105 mm**  
**Source radius: 70.5 mm**  
**Source volume: 1.640 l**  
**Source material: Water**

#### Geometry

- 1 in Stand-Off**  
1 inch Stand-off for Small containers
- 2 in Stand-Off**  
2 in Stand-Off


### Reference Configuration and Measured Efficiency

Detector  
Detector name: **GEM40-83-ICS\_55-P42871A**

Container  
Container: **20ml Scintillation Vial**

Geometry  
Geometry: **1 in Stand-Off**

Source  
Source height:  mm  
Source radius:  mm  
Source material:



## Calculated Efficiency

Calculation results

Output file: **55-P42871A 20cc Vial 1 in to 1L Marinelli.outx**

Detector name: **GEM40-83-ICS\_55-P42871A**

Container name: **1L Marinelli (132G)**

Geometry name: **No holder**

Source height: **95 mm**

Source radius: **70.5 mm**

Source volume: **1.005 l**

Source material: **Water**

Number of energies: **11**

Reference efficiency curve: **55-P42871A 20ml Vial at 1 inch**

Calculation precision: **25**

Calculation duration: **0:12**

Energy (E <sub>γ</sub> )	Effective solid angle (Ω <sub>eff</sub> )	Efficiency (ε <sub>p</sub> )
46.58	0.07859489	0.003661535
59.54	0.2328852	0.01201604
88.03	0.5012204	0.03166116
122.07	0.6333793	0.04085006
165.85	0.6875661	0.04058301
391.69	0.6875784	0.02462485
661.66	0.6659779	0.01712960
898.02	0.6452383	0.01376811
1173.24	0.6213754	0.01128761
1332.5	0.6087096	0.01026144
1836.01	0.5765179	0.008199926

Copy to Clipboard

Export to ORTEC GammaVision | Export to Canberra CAM file

Print calculation results | Print calculation report

Preview:


View efficiency curve

OK | Help

### ANGLE 4 Calculation Report


#### SUMMARY

Output file: 55-P42871A 20cc Vial 1 in to 1L Marinelli.outx  
Detector name: GEM40-83-ICS\_55-P42871A  
Container name: 1L Marinelli (132G)  
Geometry name: No holder  
Source height: 95 mm  
Source radius: 70.5 mm  
Source volume: 1.005 l  
Source material: Water  
Number of energies: 11  
Reference efficiency curve: 55-P42871A 20ml Vial at 1 inch  
Calculation precision: 25  
Calculation duration: 0:12



#### CALCULATED VALUES

Energy (E <sub>γ</sub> )	Effective solid angle (Ω <sub>eff</sub> )	Efficiency (ε <sub>p</sub> )
46.58 keV	0.07859489	0.003661535
59.54 keV	0.2328852	0.01201604
88.03 keV	0.5012204	0.03166116
122.07 keV	0.6333793	0.04085006
165.85 keV	0.6875661	0.04058301
391.69 keV	0.6875784	0.02462485
661.66 keV	0.6659779	0.01712960
898.02 keV	0.6452383	0.01376811
1173.24 keV	0.6213754	0.01128761
1332.5 keV	0.6087096	0.01026144
1836.01 keV	0.5765179	0.008199926





# ANGLE 4

## Detector Model


The Detector Model defines the physical construction of the detector. The input parameters are dependent on the detector type which may be Germanium or Sodium Iodide in Coaxial, Planar, or Well configurations. A graphic display of each model helps validate the appropriate detector type in the configuration process. Some parameters, such as the Inactive material thickness and the Contact thickness, are usually not precisely known for each detector so nominal values are typically used. These minor deviations are typically inconsequential with the Efficiency Transfer calculation method implemented in ANGLE because the minor error in transmission cancels out in the Reference and Target solid angle models. This is one of the significant advantages of Efficiency Transfer over modeling alone. And, if the detector response is affected by changes to any of these parameters, then a new Reference calibration can be generated with standard sources in the lab instead of having to return the detector to the factory for an expensive and time consuming characterization.

Detector Types:



- γ HPGe: Closed or open end coaxial, Planar, and Well
- γ Ge(Li): Closed or open end coaxial
- γ NaI: Cylinder and Well



NaI Cylinder



NaI Well



HPGe Closed End Coaxial


# ANGLE 4

## Container and Source Model

Containers define the physical holders of source or sample material, and Sources define the actual material within the container. Containers and Sources are defined independently in ANGLE to simplify the process of establishing different combinations of material and volume in each container. Common materials are pre-defined for Containers and Sources, and additional materials can easily be added based on user-defined compounds or mixtures.

Container Types:




- γ Marinelli
- γ Cylinder to define Point Source, Filter Paper, Disk, Charcoal Cartridge, and Bottles



Source	
Source height:	105 mm
Source radius:	70.5 mm
Source volume:	1.161 l
Source material:	Water

## Geometry Model

The Geometry defines the relative position of the Container to the Detector including any sample holders that may be used and up to five additional absorbing layers between the detector and the container. Common materials for sample holders and absorbers and pre-defined, and additional materials are easily added by the user.



# ANGLE 4

## Reference Efficiency Calibration

ANGLE eliminates complex, expensive, and time consuming detector characterization because the Reference Calibration can be determined by direct measurement of a known source within the lab. Optimally, the Reference Calibration is determined using a source/geometry that is similar to the one being modeled in order to minimize uncertainty in the modeled efficiency; however, any source/geometry can be used as the reference when modeling any other source/geometry with good results when all of the detector and source/geometry configuration parameters are well known.

The Reference Energy/Efficiency pairs can be manually entered into ANGLE, or imported from either ORTEC's GammaVision Efficiency Tables or Canberra's CAM files. A calibration curve is then generated using up to a 6-order logarithmic polynomial function over each of up to ten different energy regions to optimize the calibration fit. Alternatively, the reference Energy/Efficiency pairs can exclude the fit function in order to calculate the modeled efficiency for only the input energy points without any uncertainty imposed by using a fit function. The choice to use a fit function or discrete energy/efficiency pairs is largely determined by how the extrapolated efficiency calibration will be used. In many cases, the extensive calibration fit algorithms in ANGLE can achieve a much more precise calibration fit than is possible with other spectroscopy applications.

Reference efficiency curve

Experimental points  
Number of points: 10

$E_\gamma$	$\epsilon_{p,ref}$
59.54	0.0069450857
88.03	0.015166297
122.07	0.017898895
165.85	0.017524855
391.69	0.010177825
661.66	0.0071321442
898.02	0.0056748125
1173.24	0.004700603
1332.5	0.0042652902
1836.01	0.0034205941

Copy to Clipboard  
Paste from Clipboard  
Import from Canberra CAM file

Energy regions  
Number of regions: 1

[ keV	, keV ]	Polynomial orde
40	2300	6

Interpolation data

Detector  
Detector name: GEM40-83-ICS\_55-P42871A

Container  
Container: 20ml Scintillation Vial

Geometry  
Geometry: Plastic Stand-off @ 1.5 inches

Source  
Source height: 40 mm  
Source radius: 13.5 mm  
Source material: Water

Reference efficiency curve name:  
55-P42871A \_20cc Vial 1.5 in

Reference efficiency curve description:  
55-P42871A \_20cc Vial 1.5 in

New curve Load saved curve Import from ORTEC GammaVision  
OK Cancel Help

# ANGLE 4


## Calculated Efficiency

ANGLE uses an Efficiency Transfer method, which is a semi-empirical approach comprised of experimental evidence (i.e. measured efficiency of a known reference source) and mathematical comparison of effective solid angle modeling for the reference and target configurations. The precision of the effective solid angle models is based on the number of integration segments over the source volume and detector surface visible to the source, and is easily adjustable to optimize calibration accuracy versus calculation time.

The derived efficiency data can be comprised of the same energy points used in the reference calibration or user-defined energy points derived by ANGLE's robust fitting algorithms. These Energy/Efficiency pairs can then be used to generate efficiency calibrations in standard gamma spectroscopy applications. A Geometry correction file can also be generated for use in ORTEC's GammaVision application so that the final analysis results retain traceability to the Reference calibration while applying the necessary efficiency corrections to the derived geometry configuration.

Detailed and Summary reports of the reference and derived efficiency calibrations and their associated configurations are also available for verification and record retention.

$$\epsilon_p = \epsilon_{p,ref} \frac{\overline{\Omega}}{\overline{\Omega}_{ref}}$$



**Calculation results**

Output file: **55-P42871A 20cc Vial 1 in to 1L Marinelli.outx**

Detector name: **GEM40-83-ICS\_55-P42871A**

Container name: **1L Marinelli (132G)**

Geometry name: **No holder**

Source height: **95 mm**

Source radius: **70.5 mm**

Source volume: **1.005 l**

Source material: **Water**

Number of energies: **11**

Reference efficiency curve: **55-P42871A 20ml Vial at 1 inch**

Calculation precision: **25**

Calculation duration: **0:12**

Calculated values:

Energy (E <sub>γ</sub> )	Effective solid angle (Ω <sub>eff</sub> )	Efficiency (ε <sub>p</sub> )
46.58	0.07859489	0.003661535
59.54	0.2328852	0.01201604
88.03	0.5012204	0.03166116
122.07	0.6333793	0.04085006
165.85	0.6875661	0.04058301
391.69	0.6875784	0.02462485
661.66	0.6659779	0.01712960
898.02	0.6452383	0.01376811
1173.24	0.6213754	0.01128761
1332.5	0.6087096	0.01026144
1836.01	0.5765179	0.008199926

Copy to Clipboard

Export to ORTEC GammaVision | Export to Canberra CAM file

Print calculation results | Print calculation report

**ANGLE 4 Calculation Report**

**SUMMARY**

Output file: 55-P42871A 20cc Vial 1 in to 1L Marinelli.outx

Detector name: GEM40-83-ICS\_55-P42871A

Container name: 1L Marinelli (132G)

Geometry name: No holder

Source height: 95 mm

Source radius: 70.5 mm

Source volume: 1.005 l

Source material: Water

Number of energies: 11

Reference efficiency curve: 55-P42871A 20ml Vial at 1 inch

Calculation precision: 25

Calculation duration: 0:12

**CALCULATED VALUES**

Energy (E <sub>γ</sub> )	Effective solid angle (Ω <sub>eff</sub> )	Efficiency (ε <sub>p</sub> )
46.58 keV	0.07859489	0.003661535
59.54 keV	0.2328852	0.01201604
88.03 keV	0.5012204	0.03166116
122.07 keV	0.6333793	0.04085006
165.85 keV	0.6875661	0.04058301
391.69 keV	0.6875784	0.02462485
661.66 keV	0.6659779	0.01712960
898.02 keV	0.6452383	0.01376811
1173.24 keV	0.6213754	0.01128761
1332.5 keV	0.6087096	0.01026144
1836.01 keV	0.5765179	0.008199926

# ANGLE 4

## Ordering Information

Model	Description
ANGLE-BW	Advanced Gamma Spectroscopy Efficiency Calibration Software
ANGLE-UW	Update from ANGLE-B32 (Version 2 or 3) to ANGLE-BW (Version 4)
ANGLE-GW	Additional Hard Copy Documentation for ANGLE

## References

References 1 and 2 are highly recommended to the interested reader.

1. ANGLE v2.1 — New version of the computer code for semiconductor detector gamma-efficiency calculations, S. Jovanovic, A. Dlabac and N. Mihaljevic, Nuclear Instruments and Methods in Physics Research Section A, doi:10.1016/j.nima.2010.02.058.
2. Testing efficiency transfer codes for equivalence, T. Vidmar, N. Çelik, N. Cornejo Díaz, A. Dlabac, I.O.B. Ewa, J.A. Carrazana González, M. Hult1, S. Jovanovic, M.C. Lépy, N. Mihaljevic, O. Sima, F. Tzika, M. Jurado Vargas, T. Vasilopoulou and G. Vidmar, Applied Radiation and Isotopes, Volume 68, Issue 2, February 2010, Pages 355-359.
3. Reliability of two calculation codes for efficiency calibrations of HPGe detectors, K. Abbas, F. Simonelli, F. D'Alberti, M. Forte and M. F. Stroosnijder, Applied Radiation and Isotopes, Volume 56, Issue 5, May 2002, Pages 703-709.
4. Methods and software for predicting germanium detector absolute full-energy peak efficiencies, K. R. Jackman and S. R. Biegalski, Journal of Radioanalytical and Nuclear Chemistry, Volume 279, Number 1/January, 2009, Pages 355-360.
5. Calculation of the absolute peak efficiency of gamma-ray detectors for different counting geometries, L. Moens, J. De Donder, Lin Xilei, F. De Corte, A. De Wispelaere, A. Simonits and J. Hoste, Nuclear Instruments and Methods in Physics Research, Volume 187, Issues 2-3, 15 August 1981, Pages 451-472.
6. Calculation of the peak efficiency of high-purity germanium detectors, L. Moens and J. Hoste, The International Journal of Applied Radiation and Isotopes, Volume 34, Issue 8, August 1983, Pages 1085-1095.
7. ANGLE: A PC-code for semiconductor detector efficiency calculations, S. Jovanović, A. Dlabac, N. Mihaljevic and P. Vukotic, Journal of Radioanalytical and Nuclear Chemistry, Volume 218, Number 1/April, 1997, Pages 13-20.
8. On the applicability of the effective solid angle concept in activity determination of large cylindrical sources, P. Vukotic, N. Mihaljevic, S. Jovanovic, S. Dapcevic, and F. Boreli, Journal of Radioanalytical and Nuclear Chemistry, Volume 218, Number 1/April, 1997, Pages 21-26.
9. "EXTSANGLE" — An extension of the efficiency conversion program "SOLANG" to sources with a diameter larger than that of the Ge-detector, N. Mihaljevic, S. Jovanovic, F. De Corte, B. Smodiš, R. Jacimovic, G. Medin, A. De Wispelaere, P. Vukotić and P. Stegnar, Journal of Radioanalytical and Nuclear Chemistry, Volume 169, Number 1/March, 1993, Pages 209-218.
10. Introduction of Marinelli effective solid angles for correcting the calibration of NaI(Tl) field gamma-ray spectrometry in TL/OSL dating, F. De Corte, S. M. Hossain, S. Jovanovic, A. Dlabac, A. De Wispelaere, D. Vandenberghe and P. Van den Haute, Journal of Radioanalytical and Nuclear Chemistry, Volume 257, Number 3/September, 2003, Pages 551-555.
11. Contribution of  $^{210}\text{Pb}$  bremsstrahlung to the background of lead shielded gamma spectrometers, D. Mrda, I. Bikit, M. Veskovc and S. Forkapic, Nuclear Instruments and Methods in Physics Research Section A, Volume 572, Issue 2, 11 March 2007, Pages 739-744.
12. Production of X-rays by cosmic-ray muons in heavily shielded gamma-ray spectrometers, I. Bikit, D. Mrda, I. Anicin, M. Veskovc, J. Slivka, M. Krmar, N. Todorovic and S. Forkapic, Nuclear Instruments and Methods in Physics Research Section A, Volume 606, Issue 3, 21 July 2009, Pages 495-500.

Specifications subject to change  
052417