

# 828 Series

*Elemental Analysis by Combustion*



**LECO**  
EMPOWERING RESULTS

# 828 Series: Elemental Analysis by Combustion

By incorporating state-of-the-art hardware and an on-board, touch-screen software platform, the 828 Series allows you to easily handle a wide range of sample applications while significantly increasing your productivity. Available in carbon/nitrogen and nitrogen/protein configurations, the core capabilities and performance of previous generations of LECO macro combustion instruments have been maintained, while key improvements have been made in throughput, uptime, and reliability. Macro sample mass capability paired with cycle times as fast as 2.8 minutes make the 828 an ideal instrument for a diverse applications base, while delivering unparalleled sample analysis throughput.




**Maximize lab efficiency and productivity with unmatched sample throughput coupled with superior instrument uptime.**

- Rapid cycle time of 2.8 minutes
- Extended reagent lifetimes, including a reduction reagent tube lifetime of over 4,000 samples
- Rugged 30-sample position autoloader with optional expanded capacity for up to 120 samples

# User-Friendly Cornerstone® Brand Software

LECO's exclusive Cornerstone brand software with touch-screen interface enables complete access to analysis control, method settings, diagnostics, reporting, and more in a highly organized, intuitive, and immersive environment. Designed through a collaboration of customer feedback and innovative engineering, Cornerstone features all of the routine day-to-day operations within a single Analysis screen designed for speed and ease-of-use. Our innovative grouping of sample data into sets and replicates simplifies the data output and automatically calculates relevant statistics, alleviating the need for additional data processing.



Analysis

## Software Features & Benefits

The software is divided into four main sections—Analysis, Diagnostics, Settings, and Instrument—for simplified navigation and organization. Toolbars, sliders, and drop-down menus make it easy to set parameters for calibration and data processing. The software also includes real-time monitoring of ambient parameters, with fully animated system diagrams.

Advanced interactive diagnostic features include a thorough digital on-board manual, maintenance animations, photo illustrations, and screen captures that quickly provide the direction needed without having to refer to multiple manuals.


Cornerstone also supports a multilingual interface, user permissions, extended data archiving and filtering, compatibility with various Laboratory Information Management Systems (LIMS), and flexible reporting capabilities. Compliance to FDA regulations 21 CFR Part 11 for a closed analytical system is also supported.




Diagnostics > Ambients




Instrument > Furnace



Settings > Calibration




An optional Cornerstone Mobile application feature enables remote viewing of the instrument software from a smartphone, tablet, or PC. It can also be programmed to set automatic notifications from the instrument against predefined software conditions using e-mail, text message, or the Cornerstone Mobile application.



# Instrument Highlights and Features

## Application Versatility with Superior Uptime

- Quartz dual-stage furnace with exclusive oxygen environment ensures complete combustion of macro samples with maximum temperature up to 1050 °C
- Large, porous crucible facilitates macro sample combustion and extends maintenance intervals
- Thermal conductivity cell supports the flexibility of using either helium or argon as a carrier gas without a hardware change
- Dual loop aliquot doser (available in FP828 performance and CN model), provides the flexibility to optimize methods based upon sample element concentration (low/high), or analysis cost and uptime



## Low Operating Cost

- Reagent-free primary and secondary furnace eliminates the cost and downtime of furnace reagent replenishment
- Thermoelectric cooler eliminates the use of chemical desiccant reagents for the removal of combustion gas moisture
- Combustion gas aliquot system provides an extended and consistent reagent lifetime regardless of sample mass, matrix, or carbon content, including a 4000 sample reduction tube reagent lifetime

## Operator Centered Design

- Boom-mounted touch-screen user interface promotes an ergonomic workspace and optimized workflow while reducing system bench space requirements
- Open access to all reagent tubes and common maintenance areas with quick-release features speeds and simplifies preventive maintenance routines, ensuring a robust and reliable instrument with superior uptime
- *Cornerstone Mobile* remote software keeps the user updated from their smartphone on the instrument's analysis batch progress, performance, and status while away



## Reliability with Trusted Service and Support

- Knowledgeable sales force with a customer-centered focus dedicated to helping you understand and identify the best instrumentation fit for your application
- State-of-the-art Technical Service Laboratory with experienced technical application chemists to assist in method development and other application-related requests
- Global and regional LECO service network comprised of regional support centers and over 25 international LECO subsidiaries, dedicated to providing service and support offerings, including field service visits over the lifetime of the instrument

## Model Availability

Available in various models, the 828 Series is ideal for a diverse applications base and is compliant with ISO, AOAC, AACC, AOCS, and ASBC approved methods of analysis.



### FP828

The FP828 delivers fast, accurate, and precise detection of nitrogen/protein in a wide range of food, feed, and other organic matrices, with an analysis cycle time of 2.8 minutes.

- Feeds
- Petfood
- Milled products

### CN828


The CN828 delivers fast, accurate, and precise detection of carbon and nitrogen analysis in environmental and agricultural samples with an analysis cycle time of 2.8 minutes.

- Soil
- Plant tissue
- Sediments



### Organic Consumables

Get the best results from your LECO instrument by using genuine LECO consumables. Visit [www.leco.com](http://www.leco.com) for featured items, specials, and ordering information (Form number 203-828).



LECO, Cornerstone are registered trademarks of LECO Corporation.

3000 Lakeview Avenue | St. Joseph, MI 49085 | 800-292-6141 | Phone: 269-985-5496  
info@leco.com | www.leco.com | ISO-9001:2015 Q-994 | LECO is a registered trademark of LECO Corporation.

**LECO Corporation**

**LECO**

EMPOWERING RESULTS